

Agilent Stand-Alone Fluorescence Detector

Routine Fluorescence Analysis without HPLC Pump

This document describes the use of a stand-alone fluorescence detector for routine analysis without the use of the HPLC pump, TCC or column.

- 1) Remove the standard flow cell and install the standalone cell located in the fluorescence module
 - i) Unscrew the two front screws
 - ii) Remove the standard flow cell, leave all tubing connected
 - iii) Install the standalone cell and tighten the two screws
 - iv) Connect the waste tubing to the port that says "Out"

- 2) Turn on the FLD and DAD (or VWD).

- 3) Startup Chemstation and open the 'Instrument Online'

- 4) From the Method and Run Control Window
 - i) **Go to View > Method and Run Control**

- 5) Turn on detector 'Instrument > FLD > Control'

- 6) Allow thirty minutes for lamp to warm up

- 7) Setup the FLD parameters

- 8) Go to Instrument : Set Up FLD Signals
 - a) Click on 'Full'.
 - b) Click on 'Special Set points button'
 - c) Under 'Fluorescence Scan Range' box
 - (a) Set the emission range
 - (b) Set the emission range

- 9) Load sample using the syringe provided

10) Fill the sample via the 'In' valve

- i) Inject till sample is flowing through the "Out" valve
- ii) Do not remove syringe

11) Run sample (Instrument > More FLD > Take Fluorescence Scan)

12) From the Data Analysis Window

- i) Go to View > Data Analysis

13) Load signal

14) View isoabsorbance plot (Spectra > Isoabsorbance Plot)

- i) Click Exit button when finished.

15) View 3D Plot (Spectra > 3D Plot).

- i) Click Close button when finished.

16) Remove syringe and flush out cell:

- i. Initially with the sample solvent
- ii. **Finally, with methanol**

17) Flush air through the cell to remove any solvent.

18) Remove standalone cell and replace with standard flow cell

- i) Unscrew the two front screws
- ii) Remove the standalone cell
- iii) Install the standard cell and tighten the two screws
- iv) Ensure that all the tubing is correctly installed and tightly fitted

19) Turn off FLD and quit Chemstation