


Sample data

Ident Sugar_100ppb
 Sample type Sample
 Determination start 2018-10-24 10:30:23 UTC-4
 Method Sugars_Automatic Dilution
 Operator

Anions



| Peak number | Retention time min | Area (nA) x min | Height nA | Concentration mg/L | Component name |
|-------------|-----------------------|--------------------|--------------|-----------------------|----------------|
| 1 | 4.545 | 7.8303 | 70.319 | | Sorbitol |
| 2 | 5.682 | 1.1075 | 7.254 | | |
| 3 | 9.943 | 12.0899 | 55.869 | | Glucose |
| 4 | 10.832 | 10.7845 | 45.621 | | Galactose |
| 5 | 12.078 | 5.9625 | 21.096 | | Fructose |
| 6 | 16.715 | 2.5173 | 7.985 | | Lactose |
| 7 | 21.600 | 5.7315 | 10.443 | | Sucrose |